

MTN-008 Infant PK- LDMS Specimen Tracking Sheet

For login of *infant* MTN-008 stored specimens into LDMS

Participant ID		Visit Code	
<input type="text"/> <input type="text"/> <input type="text"/> - <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> - <input type="text"/> - <input type="text"/> 1	<input type="text"/> <input type="text"/> . <input type="text"/>		
Site Number	Participant Number	Chk	Who

PK SPECIMEN TIME	PRIMARY SPECIMEN TYPE	DATE COLLECTED dd-MMM-yy	TIME COLLECTED hh:mm	NUMBER OF TUBES COLLECTED	INSTRUCTIONS FOR PROCESSING LAB
------------------	-----------------------	-----------------------------	-------------------------	---------------------------	---------------------------------

INFANT BLOOD COLLECTION (PREGNANCY COHORT)					
Delivery Visit	Infant Cord Blood (CRD)			<input type="checkbox"/> Non (red top)	Process within eight hours of collection. Freeze immediately after centrifugation. Store with derivative SER
Delivery Visit	Infant Blood (BLD)			<input type="checkbox"/> Non (red top)	Process within eight hours of collection. Freeze immediately after centrifugation. Store with derivative SER

INFANT BLOOD COLLECTION (LACTATION COHORT)					
6 hour post maternal dosing	Infant Blood (BLD)			<input type="checkbox"/> Non (red top)	Process within eight hours of collection. Freeze immediately after centrifugation. Store with derivative SER

Comments: _____

Initials: _____ LDMS Data Entry Date: / / _____

Sending Staff
Receiving Staff
dd
MMM
yy
LDMS Staff

MTN-008 Infant PK- LDMS Specimen Tracking Sheet

*For login of **infant** MTN-008 stored specimens into LDMS*

Infant PK - LDMS Specimen Tracking Sheet (nonDataFax)

Purpose: This non-DataFax form is used to document collection and entry of MTN-008 infant PK blood specimens into the Laboratory Data Management System (LDMS).

General Information/Instructions: A copy of this form accompanies infant PK blood specimens (in their original specimen collection containers) to the LDMS entry laboratory. Once the specimens have been entered into LDMS, this form is kept on file at the LDMS entry laboratory. If the site chooses, a copy of this completed form may be made once the specimens have been entered into LDMS and the copy kept in the participant's study notebook. This is not required, however. Because this form is a non-DataFax form, this form should NOT be faxed to SCHARP DataFax.

Item-specific Instructions:

- **Visit Code:** Record the visit code of the visit at which the LDMS specimens were collected.
- **NUMBER OF TUBES COLLECTED:** In the box to the left of each additive type, record the total number of tubes collected. If no LDMS specimens of the primary specimen type were collected, record "0."
- **Initials – Sending Staff:** The clinic staff person who completed the form and/or who is sending the LDMS form and specimens to the LDMS entry lab, records his/her initials here.
- **Initials - Receiving Staff:** The laboratory staff person who received this form (and the LDMS specimens accompanying the form), records his/her initials here.
- **LDMS Data Entry Date:** Record the date the LDMS specimens listed on this form were entered into LDMS.
- **LDMS Data Entry Date - LDMS Staff:** The LDMS laboratory staff person who entered the specimens into LDMS, records his/her initials here.